

## Portable Digital Urine Specific Gravity Refractometer

Atago Clinical Refractometer Series

Item Number: ATA-3464



- Easy user friendly portable digital Urine Gravity refractometer
- Wide measurement range (Urine S.G. 1.0000 to 1.0600)
- Flat sensor easy to clean
- Calibration with water only
- External Light Interference (E.L.I.) function ensures accurate measurement when intense direct light is detected
- Indoor and outdoor usage
- Automatic temperature compensation
- Digital LCD

The UG-Alpha is a compact Urine Specific Gravity refractometer with digital display. Operation is simple. Just add a sample to the prism and press the START key. The measurement value of Urine Specific Gravity is displayed by using only a single drop of Urine. The measurement range has been expanded to Urine S.G. 1.0600.

### Specification

<b>Range:</b>	Urine S.G. 1.0000 to 1.0600
<b>Minimum Indication</b>	Urine S.G. 0.0001
<b>Measurement Temperature</b>	10 to 35°C (ATC)
<b>Measurement Accuracy</b>	Urine S.G. scale $\pm 0.0010$
<b>Sample Volume</b>	0.1ml
<b>Power</b>	006P Dry Battery (9V)
<b>Dimensions and Weight</b>	17 x 9 x 4 cm, 300g